



Behavior Capture

This Behavior Imaging™ application will **change the way you treat behaviors**. Capture and store behavioral events before, during, and after the event. This is more than just a video recording tool. This is a “**go-back-in-time** and see-for-myself-what-happened” application. **Highly effective and easy to use**, Behavior Capture delivers results through the ability to better understand behaviors.



Capture better data to more effectively diagnose and treat behaviors

Digital Buffering

Capture video of behavioral events in the patient’s natural environment with one click of a small remote control worn around the neck. A digital buffer will “go back in time” (up to 30 minutes) to capture the events leading up to the problematic behavior. The capture of antecedents enables more effective diagnosis and treatment of behaviors.

Annotation

Once video has been captured, it can be referenced on a computer for annotation. Users can draw attention to behaviors and further enable collaboration by noting specific instances on the video on a second by second basis. These annotations become part of the video record and can be referenced by specialists when the video is shared and played back.

Improve workflow with easy-to-use features

One-Click Capture

Small cameras are set up to blend with the natural environment. Video recording can be activated remotely with one click of a remote control worn around the neck. React quickly and capture a behavioral incident on video while it’s happening.

Auto-load to Computer

Once recording is completed, the video auto-loads to a computer desktop. No manual connection or file transfer is required to upload the video (as is the case with a standalone video or digital camera).

Portable

Behavior Capture is compact and portable for easy use in the field. There’s no need to haul around bulky video equipment that requires set up and breakdown. Use Behavior Capture for easy point and shoot video capture while you’re on the move.

Use video as data for behavior training and research

Digital File Directory

Store and manage an unlimited number of videos in a directory with an easy-to-use windows based user interface. Organize videos and create custom libraries of behaviors that can be played back for training and research.

Reporting

Behavior Capture integrates with Excel and other tools for easy reporting. Videos can be tagged with keywords to enable quick searches and grouping of videos into different categories. Behavior Connect can also be integrated with Behavior Capture for expanded record keeping and data management needs.

"Behavior Capture is the best way to do Functional Behavior Assessments...it has been a huge cost savings."

-Behavior Specialist
Cranbury, NJ

Use Cases for Behavior Capture

- Improve data collection for Functional Behavior Assessments (FBAs).
- Reduce behavior data collection errors by up to 43%.
- Diagnose and treat early developmental disorders.
- Capture and store video to document students' progress in the classroom.
- Create video glossaries for behavioral disorders.
- Capture video in home to show behaviors to specialists and doctors.
- Train staff and parents in Applied Behavior Analysis and other clinical treatment models through demonstration videos.



What's Included

- Behavior Capture Professional Software
- Logitech® 9000 USB video camera
- Radio Frequency (RF) remote control
- USB remote control receiver
- Priority phone and e-mail support
- (Optional) - pre-configured PC exceeding requirements

System Requirements

Product	Behavior Capture
Operating systems	Windows® XP or later operating system
Processor	1.5 GHz or faster processor
Memory	1 GB RAM, 20 GB of available disk space
Monitor	1024 x 768 or higher resolution
Sound cards	Windows®-compatible sound card
Ports	Minimum of two available USB ports

Maximum Potential is an authorized distributor of Behavior Imaging Solutions programs.



info@behaviorimaging.com
www.behaviorimaging.com

(c) 2010 Behavior Imaging Solutions, Inc.
Behavior Imaging is a registered trademark.
Behavior Capture and Behavior Connect are trademarks of Behavior Imaging Solutions.



MPK Contact: Garrett Butch - 877-629-7654
or Info@maximumpotentialkids.com